

Usage Guidelines for Corning® Gosselin™ Plasticware

Application Note

CORNING

Plastic Characteristics

Reagents	HDPE		PS		PP		LDPE		PET	
	68°F 20°C	122°F 50°C	68°F 20°C	122°F 50°C	68°F 20°C	122°F 50°C	68°F 20°C	122°F 50°C	68°F 20°C	122°F 50°C
Perchloric acid	G	N	G	F	G	N	G	N	F	F
Phosphoric acid 85%	R	R	R	G	R	G	R	R	F	F
Sulphuric acid 20%	R	R	R	G	R	G	R	R	N	N
Trichloroacetic acid	G	F	F	N	R	G	G	F	F	N/A
Butyl alcohol	R	R	G	G	R	R	R	R	G	G
Ethyl alcohol	R	R	F	N	R	G	R	G	G	G
Methyl alcohol	R	R	F	N	R	R	R	R	G	F
Acetaldehyde	G	F	N	N	G	N	G	N	F	N
Aniline	R	G	F	F	G	F	R	G	G	F
Benzene	G	G	N	N	G	F	G	N	F	N
Potassium dichromate	R	R	R	G	R	G	R	G	F	N/A
Bromine	F	N	N	N	N	N	N	N	N	N
Bromoform	N	N	N	N	N	N	N	N	N	N
Butadiene	F	N	N	N	N	N	N	N	N	N
Cellosolve	R	R	F	F	R	G	R	G	N/A	N/A
Chlorine (moist)	G	F	N	N	G	N	G	N	N	N
Chlorine in air	R	F	F	N	G	N	G	N	N	N
Chloroform	G	F	N	N	G	F	F	N	N	N
Amyle chloride	F	N	N	N	N	N	N	N	N	N
Ethylene chloride	G	F	N	N	F	N	G	N	N	N
Ethyl chloride	F	F	N	N	F	N	F	N	N	N
Methylene chloride	G	F	N	N	F	N	F	N	N	N
Thionil chloride	N	N	N	N	N	N	N	N	N	N
Cresol	F	N	N	N	G	F	N	N	N	N
Cyclohexane	R	G	N	N	G	F	G	F	R	N/A
P-dichlorobenzene	G	F	N	N	G	F	F	N	F	F
Dioxane	G	G	N	N	G	F	G	F	R	R
Hydrogen peroxide	R	R	R	G	R	G	R	G	R	R
Gasoline	G	G	N	N	G	F	F	N	R	G
Petroleum ether	F	N	F	F	G	F	F	F	R	G
Ethyl ether	N	N	N	N	N	N	F	N	R	G
Formaldehyde	R	R	G	G	R	G	R	G	R	R
Hexane	R	R	N	N	R	R	R	R	R	R
Mineral oil	R	R	R	R	R	R	G	N	R	R
Calcium hypochlorite	R	R	G	F	R	R	R	R	R	R
Sodium hypochlorite	R	R	R	R	R	R	R	R	G	N/A
Kerosene	G	G	N	N	G	F	F	N	R	R
Methyl isobutyl ketone	R	G	N	N	G	F	G	F	G	F
Silver nitrate	R	G	G	F	R	G	G	F	R	R
Ethylene oxide	F	F	N	N	F	F	G	F	N/A	N/A
Potassium permanganate	R	R	G	G	R	G	R	G	R	R
Phenol	G	F	N	N	G	N	G	N	N	N
Propane	F	N	N	N	N	N	N	N	F	N
Turpentine	G	G	N	N	G	F	F	N	G	G
Carbon tetrachloride	G	F	N	N	G	F	F	N	R	R
Tetrahydrofuran	G	F	N	N	G	F	F	N	N	N
Trichloroethane	F	N	N	N	N	N	N	N	N	N
Trichloroethylene	F	N	R	R	N	N	N	N	F	N
Toluene	G	G	N	N	G	F	F	N	R	N
Urea	R	R	G	G	R	R	R	R	R	R
Xylene	G	F	N	N	F	N	G	N	N	N

G = Good resistance for continuous use not exceeding 30 days

F = Fair resistance for continuous use not exceeding 7 days

N = Not recommended for continual use

N/A = Not applicable

HDPE = High-density polyethylene

LDPE = Low-density polyethylene

PET = Polyethylene terephthalate

PP = Polypropylene

PS = Polystyrene

Chemical Compatibility

G = Good resistance for continuous use not exceeding 30 days

F = Fair resistance for continuous use not exceeding 7 days

N = Not recommended for continual use

N/A = Not applicable

HDPE = High-density polyethylene

LDPE = Low-density polyethylene

PET = Polyethylene terephthalate

PP = Polypropylene

PS = Polystyrene

Reagents (20°C)	HDPE	PS	PP	LDPE	PET
Weak or diluted acids	R	R	R	R	R
Strong or concentrated acids	R	F	R	R	N
Aliphatic alcohols	R	R	R	R	F
Aldehydes	G	N	G	G	G
Bases	R	R	R	R	N
Esters	G	N	G	G	G
Aliphatic hydrocarbons	G	N	G	F	G
Aromatic hydrocarbons	G	N	F	F	N
Halogenated hydrocarbons	F	N	F	N	N
Cetane	G	N	G	G	F
Oxidizing agents	F	N	F	F	F

This data is excerpted from the Corning Gosselin Product Selection Guide (CLS-G-PSG-004) and is based on information and documentation provided by our raw material suppliers and on our own many years of technical experience. Customers should test our products to ensure that they are suitable for their own specific requirements. The information provided is for guidance only and we cannot accept further responsibility.

Warranty/Disclaimer: Unless otherwise specified, all products are for research use only. Not intended for use in diagnostic or therapeutic procedures. Corning Life Sciences makes no claims regarding the performance of these products for clinical or diagnostic applications.

For additional product or technical information, visit www.corning.com/lifesciences or call 800.492.1110. Outside the United States, call +1.978.442.2200 or contact your local Corning sales office.

Corning Incorporated Life Sciences

836 North St.
Building 300, Suite 3401
Tewksbury, MA 01876
t 800.492.1110
t 978.442.2200
f 978.442.2476
www.corning.com/lifesciences

Worldwide Support Offices

ASIA/PACIFIC

Australia/New Zealand
t 61 427286832

China
t 86 21 3338 4338
f 86 21 3338 4300

India
t 91 124 4604000
f 91 124 4604099

Japan
t 81 3-3586 1996
f 81 3-3586 1291

Korea
t 82 2-796-9500
f 82 2-796-9300

Singapore
t 65 6572-9740
f 65 6861-2913

Taiwan
t 886 2-2716-0338
f 886 2-2516-7500

EUROPE

France
t 0800 916 882
f 0800 918 636

Germany
t 0800 101 1153
f 0800 101 2427

The Netherlands
t 020 655 79 28
f 020 659 76 73

United Kingdom
t 0800 376 8660
f 0800 279 1117

**All Other European
Countries**
t 31 (0) 20 659 60 51
f 31 (0) 20 659 76 73

LATIN AMERICA
grupoLA@corning.com

Brasil
t (55-11) 3089-7400

Mexico
t (52-81) 8158-8400

CORNING | **FALCON** | **AXYGEN** | **PYREX**

For a listing of trademarks, visit www.corning.com/clstrademarks.
All other trademarks are the property of their respective owners.